


Data sheet Ingold nozzle DN 25


For installation of Jet CIP valve DN 25 on vessel.



Technical data

Operating pressure	max. 6 barg	max. 87 psi
Operating temperature	max. 150°C	max. 302 F
Material	316L	316L
Surface roughness	Ra 1.2 µm	Ra 1.2 µm
For vessel diameter	200 ... 800 mm	7.87 ... 31.5 inch

Dimensions



	D1	D2	D3	H1
[mm]	24.8	42	G 1 1/4"	46
[inch]	0.976	1.66	G 1 1/4"	1.81

(Finishing dimensions see welding instructions '0940-04 BA Ingold nozzle')

Ingold nozzle DN25	Order no.
Ingold nozzle DN 25	112568

Documentation	Order no.
Data sheet	
Welding instructions	57496

Options

Extended documentation	Order no.
Material certificate EN 10204-3.1	57497