

controllers:

Bioengineering AG 8636 Wald, Switzerland Phone +41 (0)55 256 81 11 Fax +41 (0)55 256 82 56

0946-03a RALF (Type 2017) Process specific 1/2

Bioengineering RALF – Request for quotation

Process specific configuration

Your contact information	1			
Name			X	
Company				
Designation				a constant
Address line 1				- 0 O
Address line 2				
Email				
Phone				Pier II
for Essential models instead.	h component of your RALF, pl			
Process specific configu	ıration			
Characterization of you	r process			
-		Cells	□ r	Yeast
Organism Respiration	☐ Bacteria	☐ Anaerobic	Fungi Facultative aerobic	□ feast
Other design considerations	☐ Shear sensitive	☐ Heat sensitive	☐ Highly foaming	Other
- Concrucing Consucrations	Silear SchSitive	reac sensitive	mgmy rounning	_ other
Feed				
Operation mode	Batch	Fed batch	☐ Continuous	Perfusion
Aeration				
Aeration for a microbial process	s Standard air (1 channel for 2 vvm Air with pulsed valve)	O ₂ enriched (2 channels for 2 vvm Air and 1 vvm O ₂ with pulsed valves)	☐ Anaerobic (2 channels for 2 vvm Air and 1 vvm N ₂ with pulsed valves)	Mass flow controllers Number of mass flow controllers:
Aeration for a cell process	Standard 3 gas mix	Enhanced 4 gas mix	Surface aeration	Mass flow controllers

0.1 vvm Air, 0.1 vvm O_2 and 0.1 vvm Air, 0.1 vvm O_2 ,

 0.05 vvm CO_2 , all with 0.1 vvm N_2 and 0.05 vvm



Bioengineering AG 8636 Wald, Switzerland Phone +41 (0)55 256 81 11 Fax +41 (0)55 256 82 56

0946-03a RALF (Type 2017) Process specific 2/2

Process monitoring + cor	ntrol					
Cultivation parameters Agitation, T, DO, pH incl. in standard models	OD for biomass monitoring	Exhalyzer for online gas analysis of O ₂ and CO ₂ concentration	Redox for monitoring of activity of anaerobic/micro aerobic culture	☐ Antifoam		
	☐ Balance 10-35'000 g	☐ Balance 1-16'100 g	☐ Balance 1-6'500 g	Level		
Visualization and control	■ BioSCADA (Full SCADA functionality incl. recipe configuration, data acquisition, analysis and export)	☐ Laptop (preconfigured and tested Laptop with BioSCADA)	4x 4-20 mA connected to Output: 4x 4-20 mA freely	ed to universal PID controller universal PID controller		
Dosing						
Line for	Base	☐ Acid	☐ Feed	☐ Antifoam		
	Others, please state numb	per				
Volumes + units						
Up to six units can be supplied and controlled with one PC	2.0 L: units					
	3.7 L: units					
	☐ 5.0 L: units					
	☐ 6.7 L: units					
Options						
Additional items	☐ Installation					
	☐ Installation with IQ/OQ					
	Spare parts package					
	☐ Spare filters for gas inlet and outlet ☐ Chiller 50 L, for up to 6 RALF units					
	_ clinici 30 L, for up to 6 N	TEI UIIICS				

Comments

Do you have any additional process requirements not covered by this questionnaire? Would you like a multipurpose model that can be configured for both microbial and cell culture applications? Tell us what you need and we'll put together the best version for your requirements.